



### Acrylic Brain Matrices

An economical alternative to metal brain matrices. The acrylic coronal brain matrices have the additional feature of a mid-line sagittal cut to facilitate the separation of the left and right hemispheres. All divisions are 1mm apart.

- Clean with warm water and mild detergent
- Can tolerate ethyl alcohol
- Do not autoclave

Prod #	Description	Unit
15050	Mouse, 30g Adult, Coronal, 1mm	each
15051	Mouse, 30g Adult, Sagittal, 1mm	each
15052	Rat, 175-300g Adult, Coronal, 1mm	each
15053	Rat, 175-300g Adult, Sagittal, 1mm	each
15054	Rat, 300-600g Adult, Coronal, 1mm	each
15055	Rat, 300-600g Adult, Sagittal, 1mm	each