




Mouse & Rat Stainless Steel Brain/Heart Matrices, 0.5mm

These Matrices are made of hard stainless steel. They can be cooled, cleaned, and autoclaved. They have smooth surfaces that do not chip. The investigator will be able to get repeatable coronal or sagittal sections enabling precise blocking prior to microtome sectioning.

It is advisable to not let any tissue dry onto the matrix but rather soak the matrix immediately after use for easier cleaning.

Features:

- Stainless Steel
- Highly polished
- Precision machined
- Precise blocking of the brain or heart
- Repeatable 0.5mm segments



Size

- 15045, Mouse, Coronal (40 - 75g)
- 15046, Mouse, Sagittal (40 - 75g)
- 15047, Large Rat, Coronal (300 - 600g)
- 15048, Large Rat, Sagittal (300 - 600g)
- 15056, Mouse, Heart
- 15057, Rat, Heart

A	B	C	D	E	Cavity depth
3.18	11.1	8.73	19.1	12.2	7.4
3.18	11.1	8.73	19.1	12.2	7.4
4.76	19.8	14.7	36.6	24.7	10.91
4.76	19.8	14.7	36.6	24.7	10.91
8.3	12.1				4.8
12.7	19.9				9.6

Prod #	Description	Unit
15045	Mouse, 40 - 75g, Coronal, 0.5mm	each
15046	Mouse, 40 - 75g, Sagittal, 0.5mm	each
15047	Rat, 300 - 600g, Coronal, 0.5mm	each
15048	Rat, 300 - 600g, Sagittal, 0.5mm	each
15056	Mouse Heart, 0.5mm	each
15057	Rat Heart, 0.5mm	each